

AK3L1-adenylate kinase 3-like 1 Gene

Reference: PTXBC066944

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AK3L1; AK6; AKL3L; AKL3L1; FIX; GTP:AMP phosphotransferase AK3, mitochondrial; GTP:AMP phosphotransferase, mitochondrial; adenylate kinase 3 alpha-like 1; adenylate kinase 6, adenylate kinase 3 like 1; adenylate kinase 3		

For research use only.

GENE DESCRIPTION	adenylate kinase 3-like 1
GENE ID	205
ORIGIN SPECIES	Human
ACCESSIONS	BC066944
VECTOR	pDONR223
NCBI SYMBOL	AK3L1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctccaaactcctgcgcggtcatcctcggccgccggctcgggcaaggccacgtgtgccggaggtcggcagaactttggtctccagctctcagcggccactcttgcgggagaacatcaaggccagcaccgaagtgtgagatggcaagcagtatatagaga
aaagtcttttgggtccagaccatgtgatcacacgcctaatgatgtccgagttggagaacaggcgtggccagcactggctccttgatggttttcctaggacattaggacaagccgaagccctggacaaaatctgtgaagtggatctagtgtcagtttgaatattccattgaaacacttaagatc
gtctcagccgctgttgattcaccctctagcgggaagggtatataacctggacttcaatccacctcatgtacatgtgtgatgacgtcactggtgaaccgttagtccagcaggagatgataaacccgaagcagttgctgccaggctaaagcagtcacaagacgtggcaagccagtcatt
gaattatacagagccgaggagtgctccaccaattttccggaacggagacgaacaaatctggccctacgtttacacacttttcaacaagatcacacctattcagtcacaaagaagcatattga