



Reference: PTXBC066948

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	E2-EPF; E2EPF; EPF5; ubiquitin-conjugating enzyme E2 S; E2 ubiquitin-conjugating enzyme S; E2-EPF5; ubiquitin carrier protein S; ubiquitin conjugating enzyme E2S; ubiquitin-conjugating enzyme E2-24 kD; ubiquitin-conjugating enzyme E2-24 kDa; ubiquitin-con		

For research use only.

GENE DESCRIPTION	ubiquitin-conjugating enzyme E2S
GENE ID	27338
ORIGIN SPECIES	Human
ACCESSIONS	BC066948
VECTOR	pDONR223
NCBI SYMBOL	UBE2S
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]