



Reference: PTXBC066952

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CMT1F; CMT2E; NF-L; NF68; NFL; PPP1R110; neurofilament light polypeptide; light molecular weight neurofilament protein; neurofilament protein, light chain; neurofilament subunit NF-L; neurofilament triplet L protein; neurofilament, light polypeptide 68kDa		

For research use only.

GENE DESCRIPTION	neurofilament, light polypeptide
GENE ID	4747
ORIGIN SPECIES	Human
ACCESSIONS	BC066952
VECTOR	pDONR223
NCBI SYMBOL	NEFL
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]