



Reference: PTXBC066975

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EF-hand domain-containing protein C3orf25; C3orf25; EF-hand calcium-binding domain-containing protein 12; EF-hand calcium binding domain 12		

For research use only.

GENE DESCRIPTION	chromosome 3 open reading frame 25
GENE ID	90288
ORIGIN SPECIES	Human
ACCESSIONS	BC066975
VECTOR	pDONR223
NCBI SYMBOL	C3orf25
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]