

C6orf201-chromosome 6 open reading frame 201 Gene

Reference: PTXBC066986

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C6orf201; dJ1013A10.5; protein MGC87625; chromosome 6 open reading frame 201		

For research use only.

GENE DESCRIPTION	chromosome 6 open reading frame 201
GENE ID	404220
ORIGIN SPECIES	Human
ACCESSIONS	BC066986
VECTOR	pDONR223
NCBI SYMBOL	C6orf201
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccagtggcctctgggacccggggccggccacccgctccacccctcattcccacgcctgtcacatcaagccaggagcaccggcggtagaactctgaactctctgcttagaaagaatcttagcaagacatcagttgccaaaagagattaatctgaccccaagccgaa
cagaatgccccgtggaaaagaaaaatcatcaacaatgtaactgacgggtggaaagaatgtcactgtgtgaagagaacacgaaagagcctccaatgtccaccatagttgtcagaaaactattcaaaaaatgtaccagaagacacagcctaaggaatacagcagggaactgagaaa
cctgccacaacagctaaggacacacaacagaataa