

DPY19L2P1-dpy-19-like 2 pseudogene 1 (C. elegans) Gene

Reference: PTXBC066987

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	dpy-19-like 2 pseudogene 1 (C. elegans)
GENE ID	554236
ORIGIN SPECIES	Human
ACCESSIONS	BC066987
VECTOR	pDONR223
NCBI SYMBOL	DPY19L2P1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagaacaaggagtaaacccaagccgctgcaatctccgccccagccgctgaagcggcgctacggcgctccccgcccgggagctggagggtggaagtcggccctagggcggggaaactgccagggggcgccagggaggtctccccggggaggatcccaatct
gaaaaagcgaaaaggcttgagctaaagggtggcgaagaccctctcgacccctccagttcgtcgttaattccctggcgagctccgggaagagggtcacgaactgcaggcgcggtggttccccagcagaaccactctcagcatcgccatcttggcgaattctacattggttacatt
agtaacactttttgaaaatgatcgtcatcttctcaccctcatctttggaaatgggagatgactttcgactaaatgggactttatttctactcaagaccatcattgaagcacctcattttggaggactgtggatgattgaatgacaggttactgaatatcctcttgaattaatacagtaaa
acgcttcactctttatccagaggtatcatagctgcctggtatcgacattcataggaaataatgaatttatttggactagaactaagacgtggaatgaccagaatagaacctttaatgagttcaaaagctgtgaaggattgggagatcgtgctgtttatgttggttaa