

LOC441242-hypothetical LOC441242 Gene

Reference: PTXBC066990

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized LOC441242		

For research use only.

GENE DESCRIPTION	hypothetical LOC441242
GENE ID	441242
ORIGIN SPECIES	Human
ACCESSIONS	BC066990
VECTOR	pDONR223
NCBI SYMBOL	LOC441242
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagaggccggaccagatggaattcgggggggtccgtatctctcgtcgtgagtggttgacacagcgcaacttggtccagcccgccgaggtgccaatggagatgccgatgagcttctggagctgtccgactgtcagccccgggctggccggggccctcgccac
cctcagcggcgtcgtacacccccgccaccgcccctctctctctccatcgctccgcgagtgccgcccctatgagccacagcgccgagggagccgggctcggcgctaggcaaggaaaccccgagccgggagagctggaccgggagcggccctcgcc
gctgccccttgccagggacgccaagtagagctggcagcggagctcgtccgctccgcccctcagagccgcccggcgccgggagcagaaaagccggcgccgcaaaaaagccggcgacgggggcaaaaaagccggcgccgagcgaagaagcttaggcggcgccgcaaaaaagccgc
ggcgccgagcgaagaagcttaggcggcgccgggcaaaaaagccggcgccgcaaaaaagccggcgccgcaaaaaagccggcgccgagcgaagaagcttaggcggcgccgcaaaaaagccgcggcgccgagcgaagaagcttaggcggcgccgcaaaaaagccgc
aaagccggcgccgagcgaagaagcttaggcggcgccgggcaaaaaagctgcaaaaaagccggcgccgagcgaagaagcttaggcggcgccgcaaaaaagccgcggcgccgagcgaagaagcttaggcggcgccgcaaaaaagccgcggcgccgagcgaagaagcttaggcggcgccgcaaaaaagccgc
atcacataaacatggaagtacatttgaccattttaattgtttaacattccatgtcacacgacctccacaatagaatcaatgtatataa