

ECEL1P2-endothelin converting enzyme-like 1, pseudogene 2 Gene

Reference: PTXBC067110

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			

For research use only.

GENE DESCRIPTION	endothelin converting enzyme-like 1, pseudogene 2
GENE ID	347694
ORIGIN SPECIES	Human
ACCESSIONS	BC067110
VECTOR	pDONR223
NCBI SYMBOL	ECEL1P2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcgagatcgaacgactgaaccgctgctcatgctggaggtcctcgcgagactcgggggctggacctgcgcggcgcgcgagcgcgggggttctgcgcatggacctcaactggctgcttaa