

CSDC2-cold shock domain containing C2, RNA binding Gene

Reference: PTXBC067113

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PIPPIN; dJ347H13.2; cold shock domain-containing protein C2; RNA-binding protein pippin; cold shock domain containing C2, RNA binding; cold shock domain containing C2		

For research use only.

GENE DESCRIPTION	cold shock domain containing C2, RNA binding
GENE ID	27254
ORIGIN SPECIES	Human
ACCESSIONS	BC067113
VECTOR	pDONR223
NCBI SYMBOL	CSDC2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgactcagagtcgacgtcaccctcagttgtgccccctccactcccccaagtcctcagctgtgcccacttccctccacagggagggcagcagggctctgggagcgggggtgtgtccacctcgggacctaccagccctctgcccaccaagcggaccaggacctattcagcga
cagcccggtcctcagctgccccgtgttcaaggcgctctgaagcagttctcagctcacaggccatggcttcacccccgagaaacgggtccgaggacatcttctacatgtgtctgacatcgagggggagtgacgtgccagtggaggcgacgaggtgacctacaagatgtgcct
atccctccaagaaccagaagttccagccgtggaggtgtgtctcactcagctggccccacactccccacgagacgtgtgtctgcccaggtcgtgggtcctag