



Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88
contact@proteogenix.fr
www.ProteoGenix.science

Reference: PTXBC067266

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CMT4F; periaxin		
For research use only.			

GENE DESCRIPTION	periaxin
GENE ID	57716
ORIGIN SPECIES	Human
ACCESSIONS	BC067266
VECTOR	pDONR223
NCBI SYMBOL	PRX
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1