



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

USP39-ubiquitin specific peptidase 39 Gene

Reference: PTXBC067273

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	65K; CGI-21; HSPC332; SAD1; SNRNP65; U4/U6.U5 tri-snRNP-associated protein 2; SAD1 homolog; SnRNP assembly defective 1 homolog; U4/U6.U5 tri-snRNP-associated 65 kDa protein; inactive ubiquitin-specific peptidase 39; small nuclear ribonucleoprotein 65kDa (		

For research use only.

GENE DESCRIPTION	ubiquitin specific peptidase 39
GENE ID	10713
ORIGIN SPECIES	Human
ACCESSIONS	BC067273
VECTOR	pDONR223
NCBI SYMBOL	USP39
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr