



Reference: PTXBC067299

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MDM4, p53 regulator; Mdm4 p53 binding protein homolog; MDM4-related protein 1; MDM4 protein variant Y; MDM4 protein variant G; protein Mdm4; HDMX; MDMX; MRP1; double minute 4, human homolog of; p53-binding protein; mdm2-like p53-binding protein; protein M		

For research use only.

GENE DESCRIPTION	Mdm4 p53 binding protein homolog (mouse)
GENE ID	4194
ORIGIN SPECIES	Human
ACCESSIONS	BC067299
VECTOR	pDONR223
NCBI SYMBOL	MDM4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]