

C1orf174-chromosome 1 open reading frame 174 Gene

Reference: PTXBC067302

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms:	UPF0688 protein C1orf174; RP13-531C17.2; chromosome 1 open reading frame 174		

For research use only.

GENE DESCRIPTION	chromosome 1 open reading frame 174
GENE ID	339448
ORIGIN SPECIES	Human
ACCESSIONS	BC067302
VECTOR	pDONR223
NCBI SYMBOL	C1orf174
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggagccggaagctcacaggtgcagtgccgtctcagcgcgctgaaagcacgaagttgttcggcagccaggtggcctctgccaggaagttgctgtccacgtctccaagacagcatgtctgactctctcatcccacaagccacagacagggcgaactccaagaagtcaaa
tctgacaaaaggacatcttgaagtcagaattacagaagcttgccttaagaatgacagcgtcttggccaaaagtacacctgagaccttggaaaatgagttgctgaaggcagtgccctgctccaggcagcgaagctggcgtttctgtgcagcaggggctgcaagtctctctcg
gtggctgcagagttgtgagtgactctcgttagcaaaagactagatggcctgtccgtgccaacacagtgccgggtccggagcaagaatccaacagcagctcactgtgcagaagcagaatgagccagggctacagacagagatgtgcagaagccaccactcagatggaca
acagcgtctctagatgacgacagcaatcagccaatgccctgagccggtcttggaaacgtgagctcatgcaggacctcccactgctcttcatctgtcctcaatgagcagacgagagttcagaaaaatgcattcagagccaaagatgatgatgacgacgacgatgacg
aatgtag