

C7orf49-chromosome 7 open reading frame 49 Gene

Reference: PTXBC067350

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MRI; MRI-2; modulator of retrovirus infection homolog; chromosome 7 open reading frame 49		

For research use only.

GENE DESCRIPTION	chromosome 7 open reading frame 49
GENE ID	78996
ORIGIN SPECIES	Human
ACCESSIONS	BC067350
VECTOR	pDONR223
NCBI SYMBOL	C7orf49
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaatgaggctgagatagttgatgttgcctgggaatcctgattgagagccgcaaacaggaaaaggcctgcgagcagccggccctggcggggctgataaccagagcactccctccctgctccgtgtcgcctcacacaagttctgggagcagcagtgaggaagaggacagtggg
aaacaggcactggctccaggcctcagccctccagaggccgggggttcagctctgcctgtagcaggagccctgaggaggaggaggagaggatgtgctgaaatacgtccggagatcttttcagctag