



Reference: PTXBC067435

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CPE1; CYP2E; P450-J; P450C2E; cytochrome P450 2E1; 4-nitrophenol 2-hydroxylase; CYP1IE1; cytochrome P450, family 2, subfamily E, polypeptide 1; cytochrome P450, subfamily IIE (ethanol-inducible), polypeptide 1; cytochrome P450-J; flavoprotein-linked monoo		

For research use only.

GENE DESCRIPTION	cytochrome P450, family 2, subfamily E, polypeptide 1
GENE ID	1571
ORIGIN SPECIES	Human
ACCESSIONS	BC067435
VECTOR	pDONR223
NCBI SYMBOL	CYP2E1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr