



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

NPBWR2-neuropeptides B/W receptor 2 Gene

Reference: PTXBC067482

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	neuropeptides B/W receptor type 2; G-protein coupled receptor 8; opioid-somatostatin-like receptor 8; neuropeptides B/W receptor 2		

For research use only.

GENE DESCRIPTION	neuropeptides B/W receptor 2
GENE ID	2832
ORIGIN SPECIES	Human
ACCESSIONS	BC067482
VECTOR	pDONR223
NCBI SYMBOL	NPBWR2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagcgccgctgggaccacagagcccccctgcagacagggggctctctccctccacagtggtgfcacacgtctctcaggacaatggactggccacaatgccactctccagaccactgcggtctctctatgtgctctaccgccggtgactctcggtatctgtctgtggggctga
ctggccacgcgcgcctctatctctgtaattctggggcgcccaagatgaagacgctgacacacgtctctatctgaacctggccgctgcgcagccggctctctcactgcctgactcgcctcaacatcgggagcactgcctgcagtactggccctcggggagctgctctcgaagctgtgctgtc
cgaactgcacatcaacaacatctctccagactctactctcagctgcgtagcgtggcagcgtacacttggtgtgctggccaccgtgtagtgcgcgcacatcggggcgcgaagtgcgcgcagctgtgtgtctgcgcgcgtgcctgcgcgtctacgtctctccctc
ttctctttcgtggcctgacagcaacagagctgcaggtcccaagctggggctggcggcagcaggtctcctgtcgaagcagcagccgctgtctacacgttgcctcgtggctgtgcgtccgcgtgtgcacactctgtgtctctacacagacctctgcgcagctgcggccgctg
cggctgcctcgtctgcgcaagcgtctcaggacggcgcagcggcgaagctgacacctctggctctcgtctgcgtggcgtgtgcctctcgtctgcgcaccccttcacatcgtgcctctgtgcggcctgcacacgagctgccccagaccccaactggtcatgatgtctctacgtcatcaca
gcctcagctacccaactgcacacccctctctcagctcgttctcagatgactccgaagaactccgacatattgcgcgtgca

More info: Please contact us for any question: contact@proteogenix.fr