

IL17A-interleukin 17A Gene

Reference: PTXBC067505

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CTLA-8; CTLA8; IL-17A; IL17; interleukin-17A; cytotoxic T-lymphocyte-associated antigen 8; cytotoxic T-lymphocyte-associated protein 8; interleukin 17 (cytotoxic T-lymphocyte-associated serine esterase 8); interleukin 17A		

For research use only.

GENE DESCRIPTION	interleukin 17A
GENE ID	3605
ORIGIN SPECIES	Human
ACCESSIONS	BC067505
VECTOR	pDONR223
NCBI SYMBOL	IL17A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgactcctgggaagacctcattggtgctactgctactgctgctgagcctggagccatagtgaaggcaggaatcacatccacgaatccaggatgccaaattctgaggacaagaactccccggactgtgatggtcaacctgaacatccataaccggaataccaataccaatccca
aaaggctcctcagattactacaaccgatccactcacttggaaatctccacgcgaatgaggacctgagagatataccctctgtgatctggaggcgaagtgcgcgacctgggctgcatcaacgctgatgggaacgtggactaccacatgaactctgtccccatccagcaagatcctggtc
ctgcgaggaggagcctccactgccccaaactcctcgggtggagaagatactggtgtccgtgggctgcacctgtgtcaccgccattgtccaccatgtggcctaa