

IL23A-interleukin 23, alpha subunit p19 Gene

Reference: PTXBC067511

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IL-23; IL-23A; IL23P19; P19; SGRF; interleukin-23 subunit alpha; IL-23 subunit alpha; IL-23-A; IL-23p19; JKA3 induced upon T-cell activation; interleukin 23 p19 subunit; interleukin 23, alpha subunit p19; interleukin-23 subunit p19; interleukin-six, G-CSF		

For research use only.

GENE DESCRIPTION	interleukin 23, alpha subunit p19
GENE ID	51561
ORIGIN SPECIES	Human
ACCESSIONS	BC067511
VECTOR	pDONR223
NCBI SYMBOL	IL23A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctggggagcagagctgtaatgctgtgtgctgctgcccggacagctcagggcagagctgtgctggggcagcagcctgctgactcagtgccagcagcttcacagaagctctgcacactggcctggagtgacatccactagtgggacacatggatctaagagaagagg
gagatgaagagactacaaatgatgtcccatatccagtggtgagatggctgtgaccccaaggactcagggacaacagtcagttctgcttgcgaaggatccaccagggctgatttttatgagaagctgtaggacgatatattcacagggagccttctctgctccctgatagccctgtg
ggcagcttcagctccctactggcctcagccaactcctgagcctgaggggtcaccactgggagactcagcagattccaagcctcagtcagccagcagcctggcagcgtctcctctccgcttcaaaatccttcagcagcctccagcctttgtgctgtagccgccgggtcttgcctat
ggagcagcaacctgagtcctaa