

IL4-interleukin 4 Gene

Reference: PTXBC067514

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BCGF-1; BCGF1; BSF-1; BSF1; interleukin-4; B cell growth factor 1; B_cell stimulatory factor 1; binetrakin; interleukin 4 variant 2; lymphocyte stimulatory factor 1; pitrakinra; interleukin 4		

For research use only.

GENE DESCRIPTION	interleukin 4
GENE ID	3565
ORIGIN SPECIES	Human
ACCESSIONS	BC067514
VECTOR	pDONR223
NCBI SYMBOL	IL4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggtctcactcccaactgtctccctctgtctctctgctagcatgtgccggcaactttgtccacggacacaagtgcgatatcaccttacagagatcatcaaaactttgaacagcctcacagagcagaagactctgtgcaccgagttgaccgtaacagacatcttctgctccaagaa
cacaactgagaaggaaaccttctgcagggtctgcgactgtgtctccggcagttctacagcccatgagaaggacactcgtgcctgggtgcgactgcacagcagttccacaggcacaagcagctgatccgattctgaaacggctcgacaggaacctctggggcctggcggttgtaatt
cctgtctgtgaagggaagccaaccagagtagcttggaaaactcttggaaagctaaagacgatcatgagagagaataattcaaaagtgttcgagctga