



Reference: PTXBC067536

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	G-protein coupled receptor 4; G-protein coupled receptor 19; G protein-coupled receptor 4		

For research use only.

GENE DESCRIPTION	G protein-coupled receptor 4
GENE ID	2828
ORIGIN SPECIES	Human
ACCESSIONS	BC067536
VECTOR	pDONR223
NCBI SYMBOL	GPR4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1