

ARPC3-actin related protein 2/3 complex, subunit 3, 21kDa Gene

Reference: PTXBC067747

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ARC21; p21-Arc; actin-related protein 2/3 complex subunit 3; ARP2/3 protein complex subunit p21; actin related protein 2/3 complex, subunit 3, 21kDa; arp2/3 complex 21 kDa subunit; actin related protein 2/3 complex subunit 3		

For research use only.

GENE DESCRIPTION	actin related protein 2/3 complex, subunit 3, 21kDa
GENE ID	10094
ORIGIN SPECIES	Human
ACCESSIONS	BC067747
VECTOR	pDONR223
NCBI SYMBOL	ARPC3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccggcttaccactcttctcatggatcctgatacacaactcatcggaacatggcactgttgcctatcagaagtcattcaaggacctgccccagagagacaaaagatacagatattgtgatgaagccatctattactcaaggccaatgtctcttcaaaaactatgaaattaagaatg
aagctgataggaccttgatataataactctctacatttctgaatgtctgaagaactgcaaaagtgcatttcaaaagccaagtgagaaagaatgtatagctgggaatcactaatttccattctggagagcctgttttccacttaacgaattatgccaacctgcaacaacagga
agatgaagtgatgagagcctatttacaacagctaaggcaagagactggactgagactttgtgagaagtttgcacctcagaatgataaacccagcaagtgtggacttgccttgtgaagagacagttcatgaacaagagctttcaggacctggacagtga