

RWDD4A-RWD domain containing 4A Gene

Reference: PTXBC067752

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RWDD4A; FAM28A; RWD domain-containing protein 4; RWD domain containing 4A; RWD domain-containing protein 4A; family with sequence similarity 28, member A; RWD domain containing 4		

For research use only.

GENE DESCRIPTION	RWD domain containing 4A
GENE ID	201965
ORIGIN SPECIES	Human
ACCESSIONS	BC067752
VECTOR	pDONR223
NCBI SYMBOL	RWDD4A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtccaacgaggaccaggagatggaactagaagcattacgctctatattgaaggagatgaaagtccgggaalttaagtcaggcttctttcaatagatagggtgaaatggatgccaaagccttcttaatagagattcctggacagaacatccccaaacacctccaattctat
ctatgaacgcttttttaacaacccatatacatcagctgtaaacgagagtattatgccaagctacagggaagcagtagaagctaatcttggaaccgctatgacctatacatgtttgaaatgccaaagacaataaagacagttcatggagaatcacaatcccatcaattccgcaacatcgataa
gcaatatcatctcaattgaaactcctaatacagccccatcaagtaagaaaaagacaaaaaagaacaacttcaaaagccagaagcgtaagctggcagacaaaacagatcacaaggagaacttcctcgaggctggaactgggttgatgttgtaagcatttaagcaaaactggctctaa
ggatgatgagtag