



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

TMEM183A-transmembrane protein 183A Gene

Reference: PTXBC067775

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C1orf37; transmembrane protein 183A		
For research use only.			

GENE DESCRIPTION	transmembrane protein 183A
GENE ID	92703
ORIGIN SPECIES	Human
ACCESSIONS	BC067775
VECTOR	pDONR223
NCBI SYMBOL	TMEM183A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctccggggggccgcgcctcgtacgcatcctcccccgtacatgctcccatgcccacagagagaaagcagctcaagttccggggccacgacgcctctccggcgagtgacgtgctggcgattaccaactcggatcggcgctctgtagctctgacgagtgcaaaaagcc
atagccaacgcctgttcagcagaagataaaatctcttttggccttgggaacccctcagcttctccagaggaagctcttttctggggctgtgagccctgtgacatcatcgacagcagtgatgatgatggatggatgccagagaggaaatcatcttcagagaagaactgtccagaaaaaagaaaacga
agagacacaaaagaagaactctgcacggcctctggagacagaagaagtatccatggatatttggctatcttggcctctatactctctcgtgacagattctggaattttccctgatttggcaaatccctggacactctcaactgctcactctgctccttttggcaggatgttccagagggcactacacgctg
gatctctctccctctgctctctcgcacggcgaatggaagacatcagctctctcccccctgtgctatccatgatgatgacacatttgcctgaattctcgaatctcgaagaatccctggacacttcaggaagacccccagacataaaagaattcccaaatg

More info: Please contact us for any question: contact@proteogenix.fr