



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

NMUR2-neuromedin U receptor 2 Gene

Reference: PTXBC067776

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

FM-4; FM4; NMU-R2; NMU2R; TGR-1; TGR1; neuromedin-U receptor 2; G-protein coupled receptor FM-4; G-protein coupled receptor TGR-1; growth hormone secretagogue receptor family, member 4; neuromedin U receptor 2

For research use only.

GENE DESCRIPTION	neuromedin U receptor 2
GENE ID	56923
ORIGIN SPECIES	Human
ACCESSIONS	BC067776
VECTOR	pDONR223
NCBI SYMBOL	NMUR2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr