



Reference: PTXBC067778

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DNA binding protein RFX3; transcription factor RFX3; regulatory factor X, 3 (influences HLA class II expression); regulatory factor X3		

For research use only.

GENE DESCRIPTION	regulatory factor X, 3 (influences HLA class II expression)
GENE ID	5991
ORIGIN SPECIES	Human
ACCESSIONS	BC067778
VECTOR	pDONR223
NCBI SYMBOL	RFX3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]