

MTX2-metaxin 2 Gene

Reference: PTXBC067831

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	metaxin-2; mitochondrial outer membrane import complex protein 2		
For research use only.			

GENE DESCRIPTION	metaxin 2
GENE ID	10651
ORIGIN SPECIES	Human
ACCESSIONS	BC067831
VECTOR	pDONR223
NCBI SYMBOL	MTX2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctctagtgccggaagccttcgtctccagattgcagctgcagaaccttggcctgaaatgctacattatcagcaattgaaagggagcaaaatttacttctgacaatgcagcttctcttcagctgcaggccttttgcacaaatgtgtaacttgcctatcaaagtagttgtaggcgaatgca
gaatatatgtctccatctgtaaaagtaccttttattcatgtgggaaatcaagtagtatcagaacttggccaatagtcgaattgttaaaagccaagggcattctcttagtgcagctggaggaagtcacaaagcagaatgaaagcttacatggaattagcaacaatctgtgtgactgcaga
gctgtatcttcagtggtgatgaagctacagtaggggagatcactcatgctaggtatggatctcttacccttggcctctgaatcataatttggcctatcaaaacagtgaggagtcacaaatgaagatgaaagctattggatggggaagaagactctggaccaggtcttagagatgtagacc
agtctgtcaagctctctcaaaagactgggaacacacacgtatttctcaataagcagcctactgaactgacgcactgggtatttggcctatataccattcttaccacacaaatgacaaatgatgaacttctgagaaggtgaaaaactatagcaacctccttctctgtaggagaattgaac
agcactatttgaagatcgtgtaaaaggcagcgtgcatag