



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

SRP9-signal recognition particle 9kDa Gene

Reference: PTXBC067845

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

ALURBP; signal recognition particle 9 kDa protein; signal recognition particle 9kD; signal recognition particle 9kDa; signal recognition particle 9

For research use only.

GENE DESCRIPTION	signal recognition particle 9kDa
GENE ID	6726
ORIGIN SPECIES	Human
ACCESSIONS	BC067845
VECTOR	pDONR223
NCBI SYMBOL	SRP9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccgcagtcacagacctggaggaggttcagccgcgtgccgagaagcttacctgcgtgacctatgaaggcacgtgtggtctcaaatataggcattctgatgggaactgtgtgttaaaataacagatgatttagtagacagtgctctgctctgttcctcaggctgcagtcaggtggcat
gatcatagctcactgcacctcctgacctcctgggtcctaagcgtctctctgttcctcagctcctga