

LOC388965-FUN14 domain containing 2 pseudogene Gene

Reference: PTXBC067852

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	FUN14 domain containing 2 pseudogene
GENE ID	388965
ORIGIN SPECIES	Human
ACCESSIONS	BC067852
VECTOR	pDONR223
NCBI SYMBOL	LOC388965
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccgcgtccagtcaggaaactttgagggagatattgagtcagtgacactgcagaatttctaagcagcagccgtggcgttaagctgttggccagaatctggactctcagccgaaaagtatagcgtggcaaccacactgttcattggaggtgtcactggatgggtgcacgggttt
atatccagaagggttggaaagtggctgcaacagctgtgggaggtggatttttctccttcagctgcaaacattctgggtacatcaaaagtgtggaacagagtggaagaaggacatgaagaaagccaaagacagctgaagatccctaagagcactcagatacctaaccaggtcaggag
caaagctgaggaggtgtgtcattgtgaagaaaaatgttctagtgactgggggatttttcgagggctttctgcttggcatggcatcctaa