

RAD9B-RAD9 homolog B (*S. cerevisiae*) Gene

Reference: PTXBC068031

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	cell cycle checkpoint control protein RAD9B; DNA repair exonuclease rad9 homolog B; RAD9 homolog B; RAD9 checkpoint clamp component B		

For research use only.

GENE DESCRIPTION	RAD9 homolog B (<i>S. cerevisiae</i>)
GENE ID	144715
ORIGIN SPECIES	Human
ACCESSIONS	BC068031
VECTOR	pDONR223
NCBI SYMBOL	RAD9B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgattcaaccaagattgttgcgtgatccattgttctttacatcaagtcaagaggaagtacttctgtctgtactccactgaattttgcctcaagagttcctaatgaggaatcaatggattgagcaatgctgtacacagtgagatgtttgttgctcagatgagtttgactcttcaattggaatgg
acactgagataacattttgttcaagaattgaagggaatactgacattttcagaagctacacatgctcctatatccatttttgaattccctgggaacctctgctttgagtaftgatgatagttagtggaagtaactttatttggccacattagctgatgaacaagtagagcatcttaccaca
gtcactgtgtctttcacagaacgaaaggtcagatctgattgaaaaaaggctggcaaaaatgtaactggccaggccctggaatgtatttcaaaaaagcagcaccaagaaggctttatcctaaggagactctcacaacatatctgattggaaaactgtggcagccctgcaatgaaa
gagtggatggagatgtcagtgagttatcagaagcagtgtagcaacacagagggaagtgccagggtctctgtgtctcagaagttttctgcatgttctttggagcagtttctctgaccagcaagaacactcaaccaccttcgacagcttgccaagagcaagtgacagtgaaaggaca
tgaataatgtgtgtcgcaggaagaatttaaggagtgatccaaatatttctgtattatctga