

NPM2-nucleophosmin/nucleoplasmin, 2 Gene

Reference: PTXBC068078

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	nucleoplasmin-2; nucleophosmin/nucleoplasmin 2		
For research use only.			

GENE DESCRIPTION	nucleophosmin/nucleoplasmin, 2
GENE ID	10361
ORIGIN SPECIES	Human
ACCESSIONS	BC068078
VECTOR	pDONR223
NCBI SYMBOL	NPM2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaatctcagtagcgcagtagcacggaggaaaaggcagtgacgaccgtgctctggggctgcgagctcagtcaggagagcgggactggaccttcagacccagctggaggggaagcagagctgcaggctgtgttcatacatttgcctgggggagaaagccaagaggagat
gcatcgcgtggagatctgccccagcaaaccaggaggacaagaagatgcagccgggtcaccattgcctcactccagccctcagtcctcccatggtctccatgtaggagtgagcttttccccagttactttccagctccgggctggctcaggacccgtgttcctcagtgccaggaa
cgttatgaaaaaaagctggaaaaagaagaagggaataagagccagcgttag