

NPHP3-nephronophthisis 3 (adolescent) Gene

Reference: PTXBC068082

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CFAP31; MKS7; NPH3; RHPD; RHPD1; SLSN3; nephrocystin-3; Meckel syndrome, type 7; cilia and flagella associated protein 31; nephronophthisis 3 (adolescent); nephrocystin 3		

For research use only.

GENE DESCRIPTION	nephronophthisis 3 (adolescent)
GENE ID	27031
ORIGIN SPECIES	Human
ACCESSIONS	BC068082
VECTOR	pDONR223
NCBI SYMBOL	NPHP3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggaccgctcgtcgtcgtgagccccgggcggggaagtgatcaggacacgtacggggcgggcgggcgaggcctgcgagatcccggagggtgaagccaaggcccgctgctgcgaactgttcgccgaggcgggggcgggcgagcagggcgggcgcc
cgggtcgtgcccccggggtggcgggcgggcggtgctgctggggccagcttcaagtcacgtgctcgtggtgccagagctggagtacgcggcgccgagtagcggtcaggaaggagtagagatcttgcgtcagcaagaaccaggagttgtgtccatggcgccgcg
gaggccaagctggacacggagaacaagcggctgctggggcgaactgcaggggtctcgtcgtggtggcacgatcgggtcactgcaacctggaactcctgggctcaatga