

SFT2D2-SFT2 domain containing 2 Gene

Reference: PTXBC068098

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UNQ512; dJ747L4.C1.2; vesicle transport protein SFT2B; SFT2 domain-containing protein 2; SFT2 domain containing 2		

For research use only.

GENE DESCRIPTION	SFT2 domain containing 2
GENE ID	375035
ORIGIN SPECIES	Human
ACCESSIONS	BC068098
VECTOR	pDONR223
NCBI SYMBOL	SFT2D2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacaagctgaagaaggtgctgagcgggcaggacacggaggaccggagcggcctgtccgaggtgttgaggcatcttcattaagctggagtaccagataaaaggcttcattgctgttttctataggaattctctgctcactgctgggtactgttctgtgtgggtgccaggaagg
actacacctcttcgcagtgattatacctttggtaatacgcacatcaattgggagtagcatcttcctcattgggaccagtgaaacagctgaagcgaatgttgagcctactcgtttgattgcaactatcatggtgctgtgtgttttgcaattaccctgtgttctgcttttgggtgcataacaaggacttc
acttatctctgcaatttgcagcttttggcattgacgtgtgtacagccttctcctacacattgcaaggatgctgtgaagaagtgtttgccgtgtgtcttgcataa