

TNRC6A-trinucleotide repeat containing 6A Gene

Reference: PTXBC068209

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CAGH26; GW1; TNRC6; trinucleotide repeat-containing gene 6A protein; CAG repeat protein 26; CTD-2540M10.1; EDIE; EMSY interactor protein; GW182 autoantigen; glycine-tryptophan protein of 182 kDa; trinucleotide repeat containing 6A		

For research use only.

GENE DESCRIPTION	trinucleotide repeat containing 6A
GENE ID	27327
ORIGIN SPECIES	Human
ACCESSIONS	BC068209
VECTOR	pDONR223
NCBI SYMBOL	TNRC6A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagcagggccgctgatcacattccactgaacctccctcacggaatgctctggtccgctacagttcaaaagaagagtagtgaaggcacaagaagctctgcacatgtgtgtactgggaacactactattctgtgagttgccagtgaaggagagatcagtcgtttcttgcacaaa
gccagtcttgacccttctcccggctggcagctctcgggtccagccagagccggctgggtccctcgactgttcccaactctccagccgaccgatctcaatcactggaatggtgctgggtgctgggaactaactgtggagaccttcacggcacttcactctggggaccctccat
tattccacaagcctgtgggtccccaagcagcagcgacccccgaggaattagcagcccatctcccaatcaacgcttttcttctgttgaccacctgggtgggggtggagagtcctgttaa