

LOC554223-hypothetical LOC554223 Gene

Reference: PTXBC068238

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	histocompatibility antigen-related		

For research use only.

GENE DESCRIPTION	hypothetical LOC554223
GENE ID	554223
ORIGIN SPECIES	Human
ACCESSIONS	BC068238
VECTOR	pDONR223
NCBI SYMBOL	LOC554223
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgccccgaacctgtctgtctcggggccctgtctgaccagacctgggaggcttccactcttgaggtattccacaccacatgtccggccccggcgcgcgatccccgttctctcgtggcgacgtggagacacgcagtgcggtcgacagcga
cgccacgagtcggagtgagccggaggggccgaatattgggaagaggagacagggaccgccaaggccaagcacagtgttaccgagtgaaacctgcggaccctgagcggtactacaaccagagtgaggcctgtgagtga