

PABPCP2-poly(A) binding protein, cytoplasmic, pseudogene 2 Gene

Reference: PTXBC068242

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	poly(A) binding protein, cytoplasmic, pseudogene 2
GENE ID	728773
ORIGIN SPECIES	Human
ACCESSIONS	BC068242
VECTOR	pDONR223
NCBI SYMBOL	PABPCP2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatgggatgctcctgaatggtcctcaatgtattgttgacgggttaggtctcctaaagaacagcagaacttggagctaggcgcaaaagagtcctcaagtgtaaaagtaataactgatgaagggtgaaatccaaaggatttgattgtaagctttgaaagc
atgaagatccacagacagctatggatgtgaatggaaaggagttcagtggaatcaacttacattgtctagctccgaaaaagtggagagcagacagaacttaagtcgaatttgaacagatgaagcaagatagaatccagataccagggtgctaacgtttgtgaaaaatcttgatg
agagtattgatgatgaacgtctccagaagaatttctccatttggtataatcaatagtgcaaatgtatgatggagggtgtgtgcagcaaaaggtttggctttgtatgtttctctccccagaagaagtagctaaagcagttacagcaatgaacggtaaaattgtgcccagcaagccattgcatgta
gctttagctcagtgcaagaacagcgccctggctcacttcactaacaatacatgcagaagatggcaagtgtagcagctgtgccccaccctgtaataacccctaccagccagcacctctacatgttacttcgtagcagctatcccacagaccagacacgtgctgcatactatcctctggcc
aaattgctcaactaagaccaagtcctccctgggtgttcagagtgtagacacttatccattccaatatga