



Reference: PTXBC068508

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PPP1R124; phosphatase and actin regulator 4; protein phosphatase 1, regulatory subunit 124		

For research use only.

GENE DESCRIPTION	phosphatase and actin regulator 4
GENE ID	65979
ORIGIN SPECIES	Human
ACCESSIONS	BC068508
VECTOR	pDONR223
NCBI SYMBOL	PHACTR4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1