

ENSA-endosulfine alpha Gene

Reference: PTXBC068544

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ARPP-19e; alpha-endosulfine; endosulfine-alpha variant 1; endosulfine-alpha variant 3; endosulfine alpha		

For research use only.

GENE DESCRIPTION	endosulfine alpha
GENE ID	2029
ORIGIN SPECIES	Human
ACCESSIONS	BC068544
VECTOR	pDONR223
NCBI SYMBOL	ENSA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccagaacaagaagaagaaccctgcggaggagaccggcgaggagaagcaggacacgctggagaagaaggtattctgcctgagagagctgaaggcgaagctaaaggccaataaccaagcctaggacaaaagcctggaggctccgacttctcatgaagagactc
cagaaaggcgaaggtacttgcctcaggagactacaacatggcgaagccaagatgaagaataagcagctgccaaagtcaggaccagacaagaacctggtgactggtgatcacatccccacccacagatctgccccagagaaagtctcgtctcaccagcaagcttcgggg
taa