

REEP3-receptor accessory protein 3 Gene

Reference: PTXBC068557

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C10orf74; Yip2b; receptor expression-enhancing protein 3; receptor accessory protein 3		
For research use only.			

GENE DESCRIPTION	receptor accessory protein 3
GENE ID	221035
ORIGIN SPECIES	Human
ACCESSIONS	BC068557
VECTOR	pDONR223
NCBI SYMBOL	REEP3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgtcctggtgatcctcagagccgtggtgctggtgttggatgctttatcctgcatattatcacaagctgtgaaaacaaaaacgtgaaggaaatgttcgatggatgatgtactggattgttttgcctctatactgtgattgaacagtagccgatcaaacagttgcttggttcctt
gtactatgagctgaagattgcttttgcataatggctgctttccctatacacaaggagcaagtttaatatagaaaaatccttccacttcttctcaaggaaaggagattgatgattatgtacaagcaaaaggaacgaggtatgaaaccttggttaaccttggacggcaaggtttaa
cttcagctactgctgctgttactgcagcagtaagagccaaggagcaataactgaacgtttaagaagcttcagtatgcatgtttaacaactatccaaggtgatgagcctgtgggacaaagaccataccaacctctaccagaagcaaaaaagaaaagtaaacagccccagtgatcagc
aggttatggaaatccactgaaagacggagatgagaaaacagatgaagaagcagagggccatattcagataatgagatgttaacacacaagggttcgaagatcgcaagcatgaaatctgtgaaaaccaccaaggccgcaagaggtgctgtacgggtcactaaaatacaagtg
aagaaacgaccacaagtgtatttttag