



Reference: PTXBC068558

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LERN1; LRRN6A; UNQ201; leucine-rich repeat and immunoglobulin-like domain-containing nogo receptor-interacting protein 1; leucine rich repeat neuronal 6A; leucine-rich repeat and immunoglobulin domain-containing protein 1; leucine-rich repeat neuronal pro		

For research use only.

GENE DESCRIPTION	leucine rich repeat and Ig domain containing 1
GENE ID	84894
ORIGIN SPECIES	Human
ACCESSIONS	BC068558
VECTOR	pDONR223
NCBI SYMBOL	LINGO1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr