

MGC87631-similar to hypothetical protein FLJ36492 Gene

Reference: PTXBC068597

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	coiled-coil domain containing 144 family, N-terminal like		

For research use only.

GENE DESCRIPTION	similar to hypothetical protein FLJ36492
GENE ID	339184
ORIGIN SPECIES	Human
ACCESSIONS	BC068597
VECTOR	pDONR223
NCBI SYMBOL	MGC87631
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcctcctgggtggagaaaagcggggagggctgggggtctccgaagccggcagctacgccacgaggaaagaccctagtggtggagccaggagaccagtggtacttggactaccggggaccagtggctctggctctcctacagctggggaagaacagcgtcg
gcagcgagagcaagcacggtagggcgccttagaccagctccagcacgatgcccctggagatctggcggagctccacagagctgcccggcggcggcagctccctgggtggagcacgtctggctcctggagacactggcgtggacaagaggataggagaagacatca
gcaactagttcctgaatataaggaaaaaacagacacctgaaagtctcctcaaaaataacaatccagctgctccaagccaggctgaaggagagaaggaggagtcgctgtgtacgggtggagcagatgacgtggctctgcagcttgcctcatgcagttggtgtggagatggagaccaca
gctcgactggagcgttagggagggcaccacggggcctggagagtattgcatctgcatgagcaagagttcatcaccacatcttgcagaggaaaaaagaaagggaagaatcatgtgtctaatgtgtaaggtaa