

SNX24-sorting nexin 24 Gene

Reference: PTXBC069012

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PRO1284; SBBI31; sorting nexin-24; sorting nexin 24; sorting nexin 24		
For research use only.			

GENE DESCRIPTION	sorting nexin 24
GENE ID	28966
ORIGIN SPECIES	Human
ACCESSIONS	BC069012
VECTOR	pDONR223
NCBI SYMBOL	SNX24
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggtctacatccgctcttcgctatgaagagagcgacctggagcggggatcacgggtgtttaagatagaagtgtctaatgaatggagaaaacattttgtgaaagagatacagcgaatttcattgttcacaaaaagcctaagaaatgtataaaaactccagaatcccttctaaaca
tgftaggaaactgggtcccaagcttggaaacagcgcagacaaggcttggaaacatactacaggctgtcattttagaaaatgaagaacttcccaactgtttcttgatttcctaaatgtgcgacacttgcctctctaccaaggcagaaggtgtggtatctttgatgaacagagctgaagag
tcaagcaaacgtgccaccagcctgtgctgtgttcctcaggatccatgtgtgtcctgcagccagcggtaatacaaacctgtcatctgctaacagctctttactga