

LOC285074-hypothetical protein LOC285074 Gene

Reference: PTXBC069048

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical protein LOC285074
GENE ID	285074
ORIGIN SPECIES	Human
ACCESSIONS	BC069048
VECTOR	pDONR223
NCBI SYMBOL	LOC285074
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatcagtgctgacgggtggagcgcgagctggagaaggtgctgcacaagtctcgggctacgggcagctgtgcgagcgcggcctggaggagctatcgaactacacggcggtctcaagcacgagatcctgcagagccacggccaagatgctgaattatcaggacacttcacttg
ttttgacacagtgctgtaaaagaataaaggatactgttcaaaaattggcctccgaccacaaagacatccacagcagtggtttctcgggttgaaaagccattgataaggattcactttaa