

HIST2H2BE-histone cluster 2, H2be Gene

Reference: PTXBC069193

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GL105; H2B; H2B.1; H2BFQ; H2BGL105; H2BQ; histone H2B type 2-E; H2B histone family, member Q; histone 2, H2be; histone H2B-GL105; histone H2B.q; histone cluster 2, H2be; histone cluster 2 H2B family member e		

For research use only.

GENE DESCRIPTION	histone cluster 2, H2be
GENE ID	8349
ORIGIN SPECIES	Human
ACCESSIONS	BC069193
VECTOR	pDONR223
NCBI SYMBOL	HIST2H2BE
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcctgaaccggcaaaatccgtccggccctaaaaagggtccaagaaagccgtcacaaagcccaagaagaagacggcaagaagcgcaagcgcaagagagagctactccatctacgtgtacaaggtgctgaagcaggtccacccgacacccgcatctcgtccaagg
ccatgggcatcatgaactccttctgcaacgacatcttcgagcgcgcgggagaggtctccgcctggcgactacaacaagcgtccacatcacatccgcgagatccagacggcgtgcgcctgctgctgccggcgagctggccaagcacgccgtgtccgagggcacaagg
cgggtcaccaagtacaccagctccaagtga