

DEFB4-defensin, beta 4 Gene

Reference: PTXBC069285

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DEFB4; BD-2; DEFB-2; DEFB102; DEFB2; SAP1; beta-defensin 4A; beta defensin-2; defensin, beta 2; defensin, beta 4; skin-antimicrobial peptide 1; defensin beta 4A		

For research use only.

GENE DESCRIPTION	defensin, beta 4
GENE ID	1673
ORIGIN SPECIES	Human
ACCESSIONS	BC069285
VECTOR	pDONR223
NCBI SYMBOL	DEFB4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagggtctgtatctcctctctcgttccttcattctctgatgcctctccagggtgttttgggtataggcgatcctgttacctgccttaagagtgagccatattcatccagctctttgcctagaagggtataaacaattggcacctgtggtctccctggaacaaatgctgcaaaaagcc
atga