

SCGB1C1-secretoglobin, family 1C, member 1 Gene

Reference: PTXBC069287

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RYD5; secretoglobin family 1C member 1; ligand binding protein RYD5; secretoglobin RYD5		

For research use only.

GENE DESCRIPTION	secretoglobin, family 1C, member 1
GENE ID	147199
ORIGIN SPECIES	Human
ACCESSIONS	BC069287
VECTOR	pDONR223
NCBI SYMBOL	SCGB1C1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaggggagccgtgccctcctgctgtgtggccctcacctgtttctgcatctgccgatggccacagggaggacaacgatgagttttcatggacttctgcaaacactactgggtgggacccagaggagctctatgaggggaccttgggcaagtacaatgtcaacgaagatgccaag
gcagcaatgactgaactcaagtctgagagatggcctgcagccaatgcacaaggcggagctgtcaagctgctgtgcaagtgtggcagtcaggacggtgcctaa