

WFDC9-WAP four-disulfide core domain 9 Gene

Reference: PTXBC069295

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein WFDC9; WAP9; dJ688G8.2; WAP four-disulfide core domain 9		
For research use only.			

GENE DESCRIPTION	WAP four-disulfide core domain 9
GENE ID	259240
ORIGIN SPECIES	Human
ACCESSIONS	BC069295
VECTOR	pDONR223
NCBI SYMBOL	WFDC9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagccctggattcttactcgtcatgttcacatctctggagttgtgatgcttctgctgtgctgggaagcttctggaacaaagatcccttctagatatgataagaaactgagcagtgctgggtacagccctcatataagtactgtgagaaaagggtactaaaaaatgactgtgtacgtcca
aatcatatagctgtgacactactgtggaacatctgcttagacaacgaagagcccttaaatcaatgctaaaccctag