

TSHB-thyroid stimulating hormone, beta Gene

Reference: PTXBC069298

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TSH-B; TSH-BETA; thyrotropin subunit beta; thyrotropin beta chain; thyroid stimulating hormone beta		

For research use only.

GENE DESCRIPTION	thyroid stimulating hormone, beta
GENE ID	7252
ORIGIN SPECIES	Human
ACCESSIONS	BC069298
VECTOR	pDONR223
NCBI SYMBOL	TSHB
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgactgctctcttctgatgtccatgcttttggccttgcattgtggcgaagcagatgtcttttattccaactgagtatacaatgcacatcgaaggagagagtgctctattgcctaaccatcaacaccacacatctgtgctggatattgtatgacacgggatcaatggcaactgtttctccaa
atatgctctgtccaggatgtttgacatatagagacttcactacaggactgtagaataccaggatgccacacccatgttgcctcctattttctatctgttgccttaagctgtaagtgtggcaagtgcatactgactatagtgactgcatacatgaagccatgaagacaactactgtaccaa
acctcagaagtcttattctgtaggattttctgtctaa