

RETN-resistin Gene

Reference: PTXBC069302

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ADSF; FIZZ3; RETN1; RSTN; XCP1; C/EBP-epsilon regulated myeloid-specific secreted cysteine-rich protein precursor 1; adipose tissue-specific secretory factor; c/EBP-epsilon-regulated myeloid-specific secreted cysteine-rich protein; cysteine-rich secreted		

For research use only.

GENE DESCRIPTION	resistin
GENE ID	56729
ORIGIN SPECIES	Human
ACCESSIONS	BC069302
VECTOR	pDONR223
NCBI SYMBOL	RETN
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaagctctgtctcctcctcctcctgctcctgggctgttggtgtctagcaagacctgtgctccatggaagaagccatcaatgagaggatccaggaggtcgcggctccctaataatttagggcaataagcagcattggcctggagtgccagagcgtcacctccagggggacctggc
tacttgcctcccgaggcttcgccgtcaccggctgcacttggctccgctgtggctcgtgggatgtgcgcgccgagaccacatgctactgccagtgcgcgggcatggactggaccggagcgcgctgtgtcgtgtgcagccctga