

HIST1H2AL-histone cluster 1, H2al Gene

Reference: PTXBC069306

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H2A/c; H2AFC; histone H2A type 1; H2A histone family, member C; H2A.1; histone 1, H2ai; histone cluster 1, H2ai; histone cluster 1 H2A family member i		

For research use only.

GENE DESCRIPTION	histone cluster 1, H2al
GENE ID	8332
ORIGIN SPECIES	Human
ACCESSIONS	BC069306
VECTOR	pDONR223
NCBI SYMBOL	HIST1H2AL
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgggacgcggcaagcaggaggcaagctcgcgcaaaagccaagaccgctcttctcgtgccgtctccagttcccggtggccgagtgcaaccgactgctccgcaaggcaactatgctgagcgggtcggggccggcgccggtgtacctggcggcggtgctggagtacc
tgactgccgagatcctggagctggcgggcaagcggccgcccgcgacaacaagaagaccgcaattatcccgccacttgagctggccatccgcaacgacgaggagctcaacaagctgctgggcaaaagtaaccatcgctcagggtggtgtctctgcccaacatccaggctgtgctactgc
ccaagaagaccgagagtcaccacaaggccaaggcaataa