

IL22-interleukin 22 Gene

Reference: PTXBC069308

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IL-21; IL-22; IL-D110; IL-TIF; ILTIF; TIFIL-23; TIFa; zcyto18; interleukin-22; IL-10-related T-cell-derived inducible factor; cytokine Zcyto18; interleukin 22		

For research use only.

GENE DESCRIPTION	interleukin 22
GENE ID	50616
ORIGIN SPECIES	Human
ACCESSIONS	BC069308
VECTOR	pDONR223
NCBI SYMBOL	IL22
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccgcccgcagaaatgtgagctcttctctatggggaccctggccaccagctgctcctctctcttgccctcttggtacaggaggagcagctgcgccatcagctcccaactgcaggcttgacaagtccaacttcagcagccctatatcaccaaccgcaccttcagctggctaag
gaggctagcttggtgataacaacacagacgttcgtctcattggggagaaactgtccacggagtcagtgatgagcgtgctctatctgatgaagcaggtgctgaacttcaccctgaagaagtgcgtgtccctcaatctgataggttcagccttatatgcaggaggtggtgcccttcctggc
caggctcagcaacaggctaagcacatgtcatattgaagtgatgacctgcataccagaggaatgtgcaaaagctgaaggacacagtgaaaaagcttgagagagtgaggagatcaaagcaattggagaactggattgctgtttatgtctctgagaaatgcctgcattga